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Congratulations on choosing the **eMag 34**, a pocket-size electronic magnifier designed by SCHWEIZER. SCHWEIZER's smallest electronic magnifier enables you to view printed text and images in 2 to 14 times magnification on a 3.5" display. With its simplicity, the **eMag 34** offers an enhanced design that extends operating comfort.

The **eMag 34** was designed for instant use at any location. If you have any questions or suggestions concerning the use of this product, please contact your distributor. Your feedback is greatly appreciated. We hope you enjoy working with your **eMag 34**.

Information in this document is subject to change without notice. No part of this publication may be reproduced or transmitted in any form or any means electronic or mechanical, for any purpose, without the express written permission of A. SCHWEIZER GmbH.

Copyright © 2014 A. SCHWEIZER GmbH.
All Rights Reserved.

Table of contents

What's in the box?	17
Safety and caution	18
Batterie safety notes	19
Getting to know your eMag 34	21
Getting started	22
Specifications	25

What's in the box?

The **eMag 34** package contains the following:

1. The ultra portable **eMag 34**
2. A power supply adapter
3. A wrist band
4. A user manual
5. A carrying pouch

If any items are missing from your package, please contact your local SCHWEIZER distributor.

Safety and caution

- Do not expose the **eMag 34** to extreme temperatures or humidity.
- Do not place **eMag 34** near a radiator or heater, or where it will be exposed to direct sunlight. These conditions may damage the electronics or optics of the **eMag 34**.
- Handle the **eMag 34** with care. Rough handling will damage the internal components.
- Do not remove any parts of the **eMag 34**. Contact your local SCHWEIZER distributor to service the **eMag 34** if necessary.
- Do not open any of the components as this will void the warranty. Each component contains parts which should only be serviced by an SCHWEIZER repair specialist.
- Do not use cleaning agents or abrasive materials as they will damage your unit. Avoid touching the lens of the camera and keep it clean by using a soft cloth.
- Use only the power adapter provided.
- Do not hit, vibrate, squeeze or drop the equipment.
- Do not use the unit under +5°C (41F) or above +40°C (104F).
- Do not look directly into the light source.


Batterie safety notes

We would like to inform you that our batteries/rechargeable batteries are marketed in accordance with the current law and refer you to the following points:

1. Our product contains rechargeable batteries of the type Li Ion NP60 (EE Pack).
2. To prevent harmful effects to the environment and health, used batteries/rechargeable batteries must be disposed of properly at community recycling centres or returned to the merchant. Returning batteries is required by law and is completely free-of-charge for you.
3. Batteries/rechargeable batteries should only be returned in a discharged state.
4. Valuable raw materials such as zinc, iron, and nickel can be reclaimed by recycling used batteries/rechargeable batteries.
5. When removing the battery from our devices, please observe the following procedure: turn the unit over and remove the two screws visible underneath. Remove the cover and then remove the battery.
6. The symbol of a “crossed-out trash container” indicates that it is illegal to dispose of used batteries with regular domestic waste.


The additional specification of the chemical sign for lead (Pb), cadmium (Cd), and mercury (Hg) indicates that these materials are contained in the batteries/rechargeable batteries over with the legally-specified limits.

Getting to know your eMag 34

The following figure shows the main components of the eMag 34.

1. Power on/off/ Viewing modes button
2. Power charging indicator
3. Image freeze
4. DC in
5. Centralized camera
6. Zoom button

Getting started

Power On / Off

Power On: Press and hold the green button  for 4 seconds to turn on the device.

Power Off: Press and hold the green button  again for 4 seconds to turn off the device.

After powering on the **eMag 34**, viewing mode is initially set to Full Color Mode. For optimal viewing keep the **eMag 34** approximately 2.5 cm (1") above the viewing object. Do NOT put the **eMag 34** on the object directly.

Setting the magnification

You can choose between 2 digital magnification levels by pressing the pink button . The total combined optical and digital magnification ranges from 2.0X-14X magnification.

- Users can also change the distance between the camera lens and target object for getting intermediate degrees of optical magnification.

Viewing modes

Press the green button  and select the different color modes on the LCD Screen. The **eMag 34** has three different color modes available:

- **Full Color Photo mode**

This mode will display text, pictures and objects in full color.

- **Negative reading mode**

This mode reverses the positive image. Pictures and text will be displayed as white text on a black background.

- **Positive reading mode**

This mode enhances the foreground / background contrast. Pictures and text will be displayed as black text on a white background.

Image Freeze

- Press the blue button  to freeze the image for a closer look.
- Press the button again to exit freeze image mode.
- In image freeze mode, people can still change the color modes.

Battery Indicator

The LED beside the viewing modes button is the battery indicator.

Status of the indicator:

- Red light flashing: the battery needs to be recharged;
- Green light is lit: the battery is recharging, once fully charged, the light will be turned off automatically.

Specifications

Magnification Range	2.0X - 14X
Magnification Level	2 Levels
Viewing modes	3 color modes
Image Freeze	Yes
Li-ion Battery Capacity	1500 mAh
Battery Working time	2.5 hours
Charging time	2.5-3 hours
LCD Screen size	3.5" TFT LCD
Device Dimension	105 x 72 x 17 mm 4.13 x 2.83 x 0.67 inch
Weight	120g (4Oz)
Power adapter	5V 1.5A

Manufacturer:

Aumed Group
Add. Floor 3, Building D
No. Jiada Dong Road,
Haidian Distr.
Beijing 100044
China

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