

# eMag 43 HD

Optimized image, operation and design



14-025/313118/01

# eMag 43 HD

## Technical Specifications

The portable video magnifier in optimized HD quality with folding handle can be used like a traditional magnifying glass or can also be positioned directly on the object.

- Item No. 544002
- Magnification 1,6x to 14x continuous zoom
- Viewing options Full colour (with quick change feature)  
Black on white  
White on black  
Yellow on blue  
Yellow on black  
+ 15 additional contrast modes
- Weight 228 grams
- Size 126 x 75 x 37 mm (closed)
- Operating time Approx. 2 hours
- Battery 4 rechargeable NIMH batteries AAA (approx. 2 hours charge time)
- Camera Centrally positioned HD camera
- Display 4,3 inch LCD display 16:9
- Features Snapshot function, 80 images can be saved and exported via USB  
Adjustable reading line & text mask  
LED battery indicator  
LED illumination can be deactivated  
Energy saving mode after 30 minutes
- In the box eMag 43 HD  
Case with belt loop and carrying strap  
Wrist strap  
AC power adapter  
USB cable  
Documentation
- Accessories Item 544042: Writing and reading stand

