

eMag 50 HD

The HD all-round talent
amongst portable video magnifiers



SCHWEIZER

Der LowVision-Spezialist

eMag 50 HD

Technical Specifications

The portable video magnifier in HD quality has a two-position flip-out handle and an integrated reading stand holding the screen at a natural angle.

- Item No. 550002
- Magnification 2x to 14x continuous zoom
- Viewing options Full colour (with quick change feature)
Black on white
White on black
Yellow on blue
Yellow on black
+ 15 additional contrast modes
- Weight 298 grams
- Size 139 x 89 x 20 mm (closed)
- Operating time Approx. 3 hours
- Battery Built-in rechargeable lithium-ion battery (approx. 3.5 hours charge time)
- Camera 5 megapixel HD autofocus camera
- Display 5 inch LCD display 16:9
- Features Snapshot function, 80 images can be saved and exported via USB
Adjustable reading line & text mask
LED battery indicator
LED illumination can be deactivated
Energy saving mode after 30 minutes
- In the box eMag 50 HD
Case with belt loop and carrying strap
Wrist strap
AC power adapter
USB cable
Documentation

