

Technology to enjoy...

Induction charging – how it works

The device is simply placed on the charging pad and then, as if by magic, the integrated battery is powered up without any wire connection between charging station and device. The energy is transferred by induction.



17 different magnifiers,
11 powers, 3 light temperatures:



17-016/314628-EN/00

ERGO-Lux *i* mobil – Inductive LED hand magnifier



**World
premiere!**

Wireless charging comfort!

Illuminated magnifiers with inductive charging facilitate everyday life!

The older people get, the more likely they are to be affected by various **age-related changes**: in addition to reduced visual acuity, which is often compensated by using an illuminated magnifier, **manual dexterity** also decreases gradually. Opening a battery box and replacing batteries can thus become an increasingly difficult task for older people.

SCHWEIZER has the solution: **battery changes are no longer required** to power up illuminated magnifiers.

The magic word is “induction”!*

Cut the cord – goodbye battery anxiety. The new ERGO-Lux i mobil ensures greater comfort!



Key benefits of ERGO-Lux i mobil:

- **No more changing of batteries**
The integrated battery is charged by induction.
- **At home:**
Plug in the pad once – that's it!
The magnifier has its fixed place and is always ready for use – it couldn't be easier!
- **On the go:**
Charging pad with USB port
Multiple power sources can be used (car, laptop etc.).
- **One for all, all for one...**
One pad suitable for powering up different devices.
- **3 light temperatures**
Select the colour temperature that works best for you.
- **Extra-bright lighting**
Multi-Power boost switch for 2 brightness levels.



Comfort and technology to enjoy!

*Read more how induction works on the reverse side.